Magnetic couplings, stirrer drives for chemical Reactors bmd and cyclone: Absolutely tight, safe and dependable

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
Magnetic couplings Magnetkupplungen


The main difference of a magnetic stirrer drive and a conventional stirrer drive such as rubber gland, lip seal and mechanical seal is the fact, that the internal stirrer shaft is driven without any contact to the static seal. The external magnet and the internal magnet are virtually connected through a strong magnetic field. The external magnet is located outside the pressure chamber. The internal magnet, which is inside the pressure chamber, is connected to the stirrer shaft and rotates simultaneously with the motor driven external magnet. Büchi offers low to high torque magnetic couplings and stirrer drives, which fulfill high performance mixing also in high viscous media.

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bmd magnetic coupling *bmd Magnetcupplung*

Büchi magnetic couplings hermetically seal the reaction chamber. The operation of the internal stirrer via the magnetic coupling requires a specific motor drive, which meets your requirements regarding viscosity of the media, stirrer speed and hazardous area (ATEX).

Example: Powerful stirrer drive with magnetic coupling
bmd 300, EX version, installed on polyclave.

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cyclone magnetic drive *cyclone Magnetantrieb*

Büchi Magnetic drives incorporate a magnetic coupling with an electric motor. The stirrer speed is controlled by its own controller unit.

Example: Powerful cyclone 300 magnetic stirrer drive with stirrer servo motor controller.

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bmd 800


bmd 1200

For technical Details of standard bmd 075 to bmd 5400 please refer to the table on reverse side

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Büchi cyclone 075 as well as the cyclone 300 are complete stirrer drives. They contain a magnetic coupling bmd 075 / bmd 300 with an electric motor. They can be installed on the reactor cover plate with an M36x1.5 thread.

The cyclone 300 controller enables stepless stirrer speed adjustment of the servo motor. Stirrer speed and torque can be transferred via a RS 232 interface to bds measuring / display unit and/or to the bls software. The stirrer speed can be set electrically by a 4...20 mA Signal.

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cyclone 250 «i», bmd 250 «i»

The inert magnetic coupling bmd 250 «i» and magnetic stirrer drive cyclone 250 «i» are specially developed to avoid corrosion in processes with highly aggressive media, or for processes where contact with metal parts is not permissible. Developed on the proven Büchi magnetic coupling design, non-metallic materials such as ceramic and PEEK are used for all components in contact with the process media. High mechanical and chemical stability guarantees trouble-free operation under the harshest conditions. This new technology allows us to provide complete non-metallic pressure reactors with powerful high-torque and leak-free agitators.

cyclone «i»

The cyclone «i» is an economic inert stirrer concept with capability of mixing low viscosity media only. It can be used at pressures up to 6 bar in the inertclave.

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